

## Verification of Proficiency Module for Chest Tube Insertion SIU Surgical Skills Laboratory

*OBJECTIVES: By the end of this module the resident will meet proficiency standards for the performance of:*

- Chest Tube Insertion

### GUIDELINES FOR PRACTICE

- View the videotape of the expert performance of a chest tube insertion
- Using the performance check list as a guide, practice on the SIU chest wall model focusing on efficiency and economy of motion by observing the time it takes you to complete the procedure

### INSTRUCTIONS FOR VERIFICATION FOR PROFICIENCY TESTING

When you believe your practice sessions have resulted in proficiency in this skill area, you may request from the skills lab personnel to be videotaped for your verification of proficiency performance. Once this videotape has been completed, it will be submitted to a faculty member for evaluation. Their decision will either verify your proficiency in this skill or recommend further practice.